

What should I do if my data is not accurate?

1. On healthmov homepage, tap Settings on the top righthand corner
2. Go to "Connections."
3. Ensure that your device is successfully connected.
4. Open your device app to sync your data from your device to the app.
5. Return to healthmov to refresh the data.

If you still notice a discrepancy, please email us at support@healthmov.com. Include screenshots from both your device app and healthmov app showing the discrepancy for specific days.

Revision #4

Created 6 June 2024 12:04:42 by Saloni

Updated 9 September 2024 13:22:03 by Saloni